

IKA

designed for scientists

HRC 2 GREEN control



Operating instructions

EN

Device setup




Fig. 1

Item	Designation
1	Wireless Controller (WiCo)
2	Handle
3	Venting grid
4	Filling opening lid
5	RS 232 port
6	USB port
7	External temperature sensor socket
8	Multifunction port
9	Overflow
10	Power socket
11	Power switch
12	Pump connection OUT
13	Pump connection IN

Contents

	Page
Device setup	2
EU Declaration of conformity	4
Note for USA (FCC)	4
Note for Canada (IC)	4
Explication of warning symbols	4
Safety instructions	5
General information	5
Disposal of device.....	6
Fluids.....	6
Battery pack RB 1 (for WiCo).....	6
Intended use	7
Use	7
Range of use	7
Wireless remote control.....	7
Unpacking	8
Unpacking	8
Delivery scope	8
Preparations	8
Setting up	8
Connecting the tubings	8
Inserting battery into the WiCo	9
Mounting the WiCo to the station	9
Filling and draining	9
Fluid (Standard information for IKA fluid)	11
Moving the device.....	12
Charging the Battery Pack RB 1 (rechargeable battery)	12
Changing battery to WiCo	12
WiCo – holder WH 10	12
Operator panel and display	13
station	13
WiCo	13
Setting the safety temperature	14
Commissioning	14
Useful information	15
Working with WiCo	16
Working screen at the time of delivery.....	16
Explanation of symbols on the working screen.....	16
Menu navigation and structure	17
Menu (Details).....	19
Interface and output	25
Maintenance and cleaning	29
Error codes	30
Accessories	31
Technical data	32
Warranty	34
Pump performance curve	34

EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the regulations 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the standards or other normative documents: EN 61010-1, EN 61010-2-010, EN 61010-2-011, EN 61326-1, EN 60529, EN ISO 12100 and DIN 12876-1, -2, -3.

Bluetooth® module:

Directive: 2014/53/EU

Standards: EN 300328, EN 301489-17, EN 301489-1, EN 60950-1

A copy of the complete EU Declaration of Conformity or further declarations of conformity can be requested at sales@ika.com.

Note for USA (FCC)

This equipment complies with Part 15 of the FCC rules. Any changes or modifications not expressly approved by the Manufacturer could void the user's authority to operate the equipment. This device complies with Part 15 of the FCC rules subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept all interference received, including interference that may cause undesired operation.

NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Note for Canada (IC)

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

This device complies with Health Canada's Safety Code 6/IC RSS-210. The installer of this device should ensure that RF radiation is not emitted in excess of the Health Canada's requirement.

Explication of warning symbols



Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.



Indicates a hazardous situation, which, if not avoided, can result in death, serious injury.



Indicates a potentially hazardous situation, which, if not avoided, can result in injury.



Indicates practices which, if not avoided, can result in equipment damage.



Indicates the risk of fire or explosion.



Indicates crushing risk to fingers/hand.



Indicates a hazardous situation that cause by a hot surface.



Warning: Flammable material.



Warning: Hot surface.

Safety instructions

General information:

- **Read the operating instructions completely before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where it can be accessed by everyone.
- Ensure that only trained staff work with the device.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- The device must only be used in a technically perfect condition.

NOTICE

- Pay attention to the marked sites in Fig. 1.

Device setup:

CAUTION

- The power switch of the device must be accessed immediately, directly and without risk at any time. If access to the power switch cannot be ensured, an additional emergency stop switch that can be easily accessed must be installed in the work area.
- Insufficient ventilation may result in the formation of explosive mixtures. Therefore, only use the device in well-ventilated areas and keep all ventilation openings open.

WARNING

- Beware of the high dead weight of the device when transporting.
- Ensure that your fingers do not get crushed when setting down the device.
- Because of the heavy weight of the device, at least two person are needed for carrying the device.

NOTICE

- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.

Working with the device:

DANGER

- Do not use the device in explosive atmospheres, it is not EX-protected.
- With substances capable of forming an explosive mixture, appropriate safety measures must be applied, e.g. working under a fume hood.
- To avoid body injury and property damage, observe the relevant safety and accident prevention measures when processing hazardous materials.

WARNING

- Process pathogenic material only in closed vessels under a suitable fume hood. Please contact IKA application support if you have any question.
- When using critical or hazardous materials in your processes, IKA recommends to use additional appropriate measures to ensure safety in the experiment. For example, users can implement comprehensive monitoring equipment.
- The safety circuit (safe temperature) must be adjusted so that the maximum permissible temperature cannot be exceeded even in the event of a fault. Check the safe temperature circuit on a regular basis (see section "Setting the safety temperature").

CAUTION

- Do not start up the device if:
 - It is damaged or leaking
 - Cable (not only the power cable) is damaged.
- Power cable should not get in contact with hot parts and fluids.
- Be careful when filling a hot bath.

CAUTION

- At high operating temperature, the temperature of housing parts, surfaces and tubes can exceed 70 °C.

NOTICE

- Covers or parts that can be removed from the device without tools must later be refitted to ensure safe operation. This will prevent the infiltration of foreign objects, liquids and other contaminants.
- The feet of the device must be clean and undamaged.
- Always empty the bath if the unit is unused for long time.
- A laboratory circulator heats/refrigerates and circulates fluid according to specified parameters. This involves hazards due to high and low temperatures and general hazards due to the device of electrical energy. The user safety can not be ensured simply with design requirements on the part of the device. Further hazard sources may arise due to the type of tempering fluid, e.g. by exceeding or undercutting certain temperature thresholds or by the breakage of the container and reaction with the carrier fluid. It is not possible to consider all eventualities. They remain largely subject to the judgment and responsibility of the operator. For this reason, it become necessary for user to take other precautionary safety measures.
- The device may only be used as prescribed and as described in these operating instructions. This includes operation by instructed specialist personnel.

- When device is used for external circulation, extra precaution must be taken for hot / cold fluid leakage due to damaged hose.
 - Use suitable hoses for connection.
 - Secure hoses and tubes against slippage and avoid kinks.
 - Check hoses, tubes and bath at regular intervals for possible material fatigue (cracks/leaks).
- Do not transport or empty the bath while it is still hot/cold. This may result in accidents, especially scalding injuries.
- Always empty the bath before moving the unit.

Accessories:

- Protect the device and accessories from bumping and impacting.
- Check the device and accessories beforehand for damage each time when you use them. Do not use damaged components.
- Safe operation is only guaranteed with the accessories described in the "Accessories" section.

Power supply / switching off the device:

- The voltage stated on the type plate must correspond to the power voltage.
- The device can only be disconnected from the power supply by pulling out the power plug or the connector plug.
- The device must only be operated with the original power cord set.
- The socket for the power cord must be easily accessible.
- Socket must be earthed (protective ground contact).
- After a power failure during operation, the device may start automatically (depending on operating mode).
- Disconnect the power plug before attaching or changing any accessories.
- Disconnect the power plug before cleaning, maintenance and transportation of the unit.

Maintenance:

- The device must only be opened by trained specialists, even during repair. The device must be unplugged from the power supply before opening. Live parts inside the device may still be live for some time after unplugging from the power supply.

Disposal instructions:

- The device, accessories and packaging must be disposed of in accordance with local and national regulations.

Fluid:



WARNING

- Only use the fluids, which fulfill the requirements for safety, health and device compatibility. Be aware of the chemical hazards that may be associated with the bath fluid used. Observe all safety warning for the fluids.

- Depending on the bath fluid used and the type of operation, toxic can arise. Ensure suitable extraction.
- Do not use any fluid that may cause hazardous reactions during processing that pose significant risks to health, safety and the environment.
- The safety temperature must be set in accordance with EN 61010-2-010 chapter "Requirements for equipment containing or using flammable liquids".
 - The surface temperature of the flammable medium that is exposed to air should not exceed its flash point. A danger usually arises if a medium is heated in open vessels.
 - The surface temperature of the heating element (e.g. the heating plate of a magnetic stirrer and/or heater of a thermostat) at the surface of the flammable medium and in contact with air should not exceed the value of $(t - 25) ^\circ\text{C}$ (= set value of the safety circuit), whereby **t** is the fire point of the medium. In case of uncertainties regarding the fire point, we recommend using the lower flash point as **t**.
A danger (glass breakage) usually arises if a medium is heated in glass vessels .
- If a setting made by the user (medium temperature or safety temperature) could bring a flammable medium into a state in which the conditions mentioned above could be exceeded, additional measures must be introduced that will protect the user from this danger.
- Examine regularly the function of the safety temperature limiter.



CAUTION

- Beware of the risk of burning due to delay in boiling!



NOTICE

- Only use recommended bath fluid. Only use non-acid and non corroding fluid.
- Never operate the device without sufficient fluid! Check the fluid level detection at a regular basis.
- If water is used at higher temperature, there is heavy loss of fluid due to the evaporation of the steam. When changing the bath fluid type, must remove the remaining water from the complete system (including hoses and external devices). When doing this, also open the stopper and union nuts caps of the pump outputs and inputs and blow compressed air through all the pump outputs and inputs.
- Continuous monitoring of the filling level of the bath fluid is required, especially at high temperatures.
- To ensure a sufficient fluid circulation, the viscosity of the bath fluid must not exceed of 50 mm²/s at the lowest operating temperature.
- Untreated tap water is not recommended. It is recommended to use distilled water or high purity water (ion exchangers) and add 0.1 g soda (sodium carbonate Na₂CO₃) /liter, to reduce corrosive properties.

- Don't use following fluids:
 - Untreated tap water
 - Acids or bases
 - Solutions with halides: chlorides, fluorides, bromides, iodides or sulfur
 - Bleach (Sodium Hypochlorite)
 - Solution with chromates or chromium salts
 - Glycerine
 - Ferrous water.

Refrigerant:



- Flammable refrigerant used:
 - Do not use mechanical devices or other means to accelerate the defrosting process, other than those recommended by the manufacturer.
 - Do not damage the refrigerant circuit.
 - Only the manufacturer may open and repair the closed refrigeration system.
- The minimum room size is 10 m³, so that in the event of an error (exit of the complete refrigerant filling quantity) no flammable refrigerant-air mixture can arise (R 290 max. 8 g/m³ ambient air).
- The device has been constructed to the requirements of EU (European Union) and EFTA (European Free Trade Association) countries.

Battery pack RB 1 (for WiCo):

- If during operation the battery pack RB 1 (rechargeable battery) becomes fully discharged, the device (**station**) will continue to run or is shut down depending on the value settings for "Time Out", "Safe Speed" and "Safe Temperature". If the device is set so that it continues to run when the battery of the **WiCo** is fully discharged, the only means of switching the **station** off are the "safe STOP", "ON/OFF" key and the "Power switch"!



- Please note the following safety instructions for the battery pack RB 1 (rechargeable battery):
- Keep the battery pack out of reach of children at all times.
- Store the battery pack in a cool, dry place.
- Never throw the battery pack into a fire. Keep it away from direct sunlight and temperatures above 60 °C. High temperatures will damage the battery pack and render it unusable. Temperatures above 100 °C may cause it to explode.
- Never throw the battery pack into water or expose it to moisture. Water may lead to a short-circuit, causing the battery pack to explode.
- Do not deform or crush the battery pack or damage it in any other way. This can cause battery fluid to leak and/or the battery pack to explode.

- When not in use, keep battery packs away from paperclips, coins, keys, nails, screws or other small metal objects which could cause the contacts to be bridged. Short-circuiting may result in an explosion.
- Explosion of a battery pack may release battery fluid and cause a fire.
- The lithium polymer battery pack must only be used and charged in IKA products designed for use with this battery pack.
- When the battery pack is inserted it should slide in easily and without resistance. Do not force it.
- If the battery pack is removed for an extended period of time, store it in a sealed plastic bag to prevent short-circuiting due to moisture or contact with metal components.
- The operating temperature range of the battery pack is from 0 °C to + 45 °C. Note that the battery pack capacity will be reduced at temperatures below 20 °C.
- Only the rechargeable battery types recommended in the technical data may be used in the device!



Do not charge batteries that have leaked or that are discolored, deformed or damaged in any other way.

Disposal instructions:

- When disposing of the IKA battery pack, please tape over the contacts with adhesive tape to prevent short-circuiting due to moisture or contact with metal components. Short-circuiting may result in an explosion.
- Do not throw used battery packs into your household waste. Dispose of them properly in accordance with statutory regulations.



End users are obliged by law to return all used disposable and rechargeable batteries. Throwing them into the household waste is prohibited. Disposable/rechargeable batteries containing harmful substances are marked with this symbol to indicate that they may not be disposed of as household waste.

- You can return used disposable and rechargeable batteries free of charge to your local authority collection site or to any battery retailer. In doing so you will be complying with statutory regulations and helping to protect the environment.
- Batteries must be disposed of in accordance with local and national regulations.

Intended use

Use:

Use **HRC 2** (Heated Recirculating Chiller) for heating/refrigerating and circulating fluids.

Intended Use: Tabletop device

Range of use:

Indoor environments similar to that a laboratory of research, teaching, trade or industry.

Wireless remote control:

Before using the wireless link between the **WiCo** and the laboratory device, first check whether your region is included in the radio communication approval for the device. If it is not, remote control can also be performed using a USB cable.

The safety of the user cannot be guaranteed:

- If the device is operated with accessories that are not supplied or recommended by the IKA.
- If the device is operated improperly or in contrary to the IKA specifications.
- If the device or the printed circuit board are modified by third parties.

Unpacking

Unpacking:

- Unpack the device carefully.
- Any damage should be notified immediately to the shipping agent (post office, railway network or logistics company).

Delivery scope:

- **HRC 2 GREEN control with WiCo**
- Power cables
- Hose olive DN 8 (2 pieces) see Fig. 2
- Hose olive DN 12 (2 pieces) see Fig. 2
- Screwdriver (use for safety circuit) see Fig. 3
- **WiCo holder WH 10**
- **OS 1.0** power supply (for **WiCo**) see Fig. 4
- Temperature sensor **Pt 100.30**
- USB 2.0 cable micro A–micro B
- USB 2.0 cable A–micro B
- Plastic cap (for “**Overflow**” connector)
- User guide
- Warranty card.




Fig. 2




Fig. 3




Fig. 4

Preparations

Setting up:

- Place the unit on an even, stable, clean, nonslip, dry and fireproof surface.
- Keep at least 20 cm of open space on the front and rear side.
- The place for installation should be large enough and provide sufficient air ventilation to ensure the room does not warm up excessively because of the heat from device radiates to the environment.
- Do not set up the device in the immediate vicinity of heat sources and do not expose to sun light.
- Cooling machine, pump motor and electronics produce intrinsic heat that is dissipated via the venting grids (3)! Never cover these venting grids!

Note: After setting up the device, wait at least one hour before starting the operation to avoid the damage to the cooling system.

Connecting the tubings:

Unscrew the union nuts and stoppers using a wrench (AF 19) from the pump connector **IN** (13) and **OUT** (12).




Fig. 5

Connect the hoses for circulating the external system to the pump connection M 16 x 1 for **IN** (13) and **OUT** (12) directly or with the olives.

Screw the hose olives to the pump connection **IN** and **OUT** with union nuts. Slide the hoses (DN 12) onto the olives. The hoses must be secured with suitable clamps.




Fig. 6

Via a hose to the **"Overflow"** connector, overflowing fluid can be directed into a suitable vessel. The vessel should be positioned lower than the **"Overflow"** connector.




Fig. 7

Note: Please note the permissible temperature range of hoses. For hot fluids we recommend the IKA **LT 5.30/LT 5.31/ LT 5.32** hoses.

When the external system is not necessary, please seal the pump connection **IN** and **OUT** with the existing union nuts and stoppers.

In case of closed external system, if there is in any case no risk of exceeding the maximum fluid level, you can close the **"Overflow"** connector with the included overflow cap. This minimizes the energy loss and protects, in case of long-term applications with low temperatures and high humidity, against unwanted water input and the associated volume increase.

Check fluid heat expansion!

Inserting battery into the WiCo:





Fig. 8


Mounting the WiCo to the station:



Note: If the **WiCo** need be permanently attached to the **station**, we strictly recommend to fasten the unlocking button with the integrated screw (turn counter clockwise).

Filling and draining:

Before filling the fluid into the bath, open the front cover as indicated in following figure.



Check and make sure that the drain valve is closed (rotate clockwise to the stop position, see Fig. 11).


Note: Please note information in section "Commissioning".

Connect the power plug and turn on the device with power switch (11).

The low level icon () appears on the display of the **WiCo**. Meanwhile, the bottom LED segment of fluid level indicator lights up in red color to indicate the low fluid level.

Open the filling opening lid (4) and remove the stopper from the filling opening. Then, fill fluid to the bath.

Note: The stopper must be kept being plugged in during during operation except for filling fluids.



Note: Pay attention to the fluid level information.

Fluid level information on the **WiCo** display:

- Low level
- High level

Fluid level information on the fluid level indicator:




Fig. 12

To drain the fluid from the bath, connect a hose to the drain port and turn the drain valve in counter clockwise direction with a straight screwdriver.




Fig. 13

Warning: Do not empty the bath while it is still hot/cold, there is a risk of burning/freezing.


Fluids (Standard information for IKA fluid):

No.	IKA Designation	Operating temperature range for open bath application (°C)	Operating temperature range for closed bath applications (°C)	Safety temperature (°C)	Flash point (°C)
0	CF.EG28.N10.80.8	-10 ... 80	-10 ... 80	90	115
1	CF.EG39.N20.80.16	-20 ... 80	-20 ... 80	90	115
2	CF.EG44.N25.80.19	-25 ... 80	-25 ... 80	90	115
3	CF.EG48.N30.80.22	-30 ... 80	-30 ... 80	90	115
4	UF.Si.N30.150.10LV	-30 ... 130	-30 ... 150	145 ❶	> 170
5	HF.Si.20.200.50	20 ... 200	20 ... 200	255	> 280
6	HF.Si.20.250.50A	20 ... 200	20 ... 250	255	> 280
7	Water ❷	5 ... 95	5 ... 95	-	-
8	Customized ❸				

The listed fluid temperature range may deviate from the device operating temperature ranges and settings.

Nomenclature for IKA fluids:

CF.EG28.N10.80.8 --



(1) Classification:

HF: Heating Fluid
 CF: Cooling Fluid
 UF: Universal Fluid

(2) Chemical composition:

Si: Silicone oil
 EG: Ethylene Glycol

(3) Temperature range: (Minimum temperature. Maximum temperature)

N: Negative Temperature

(4) Viscosity:

Viscosity at 25 °C for Heating Fluid (HF)
 Viscosity at -20 °C for Cooling Fluid (CF)
 Viscosity at 25 °C for Universal Fluid (UF)

Dynamic viscosity [mPa•s] is a product of kinematic viscosity [mm²/s] and density [kg/m³] of the fluid divided by 1000.

(5) Additional information:

A: Oil Additives
 LV: Low Viscosity

❶ **Note:** for open bath application!

❷ **Note:** Tap water may be unsuitable for operation because the calcium carbonate content may cause calcification. High purity water (from ion exchangers) and distilled or bi-distilled water are unsuitable for operation due to corrosive properties of these media. High purity water and distillates are suitable as a medium after adding 0.1 g soda (Na₂CO₃, sodium carbonate) per liter of water.

❸ **Note:** The temperature limit values are adjustable in accordance with the fluid used.

Moving the device:

Empty all fluid in the bath before moving device from one place to other place.

The device can be lifted up and moved by holding the top handles or bottom handles.

It can also be moved on flat surface by lifting and pushing the front of the device. It is easy to move the device with the help of the wheels.

The angle of inclination should never be more than 10 ° at any direction when move the device!

Note: The device must not be moved during operation. After moving the device, you must wait at least one hour before restarting the unit.




Fig. 14

Charging the Battery Pack RB 1 (rechargeable battery):

The battery pack in the **WiCo** can be charged by any of the following means:

- On the **station**
- Via a USB cable at the PC or **station**
- Via an **OS 1.0** power supply unit.

Changing battery to WiCo:




Fig. 15

Note: Please observe the relevant safety instructions in the "Safety instructions" section for the RB 1 battery pack!

WiCo – holder WH 10:




Fig. 16

Operator panel and display

station:




Fig. 17

Item	Designation	Function
A	On/Off button:	Switch on/off the station .
B	Adjustable safety circuit:	Adjust the safety temperature limit with delivered screwdriver.
C	LED bar:	Display different status of the circulator with different color.
D	Contacts:	Communicate and charge the WiCo .
E	"safe STOP" button:	Safe stop the working of the station in emergency.
F	Bluetooth® searching button:	Search the WiCo when the Bluetooth® is active.
G	Bluetooth® LED:	Bluetooth® indicator.
H	Fluid level indicator:	Indicate the bath fluid level (when the fluid level is too low or too high, the bottom or the top LED segment change into red color).

WiCo:



Fig. 18

Item	Designation	Function
L	On/Off button:	Switch the WiCo on and off.
M	"int/ext" button:	Switch between the internal and external temperature display and control.
N	"Timer/Pump" button:	Switch between the timer and pump display.
O	Graph button:	Display time/temperature graph.
P	Rotating/pressing knob:	Set the temperature value.
R	Rotating/pressing knob:	Start/Stop the heating/cooling function. Set the pump speed value. Start/Stop the pump function.
S	Display:	Screen.
T	"Back" button:	Return to the previous menu level.
U	"Menu" button:	Press it once: main menu is displayed. Press it a second time: back to the working screen.
V	Lock button:	Lock/Unlock the knobs and buttons.

Setting the safety temperature

Setting the safety temperature with screwdriver delivered with the device.




Fig. 19

The safety temperature setting will appear on the display.

Factory setting: maximum value.


Adjustment range: 0 to 110 °C.

Warning: The safety temperature must always be set to at least 25 °C lower than the fire point of the fluid used. In case of uncertainties regarding the fire point, we recommend setting the safety temperature to at least 25 °C lower than the flash point of the fluid used.

Commissioning

Note: Before commissioning, make sure that the device has not been moved in one hour!

Check whether the voltage given on the type plate corresponds to the available power voltage.



The socket used must be earthed (fitted with earth contact).

If these conditions have been met, the machine is ready for operation when the power plug is plugged in.

If these conditions are not met, safe operation is not guaranteed and the machine could be damaged.

Observe the ambient conditions (temperature, humidity, etc.) listed under "Technical Data".

After pressing the power switch (11) of the **station**, the screen of the **WiCo** displays the device designation and the software version after a beep.



Fig. 20

After several seconds, screen display the information of **WiCo**.

Information	
Internal Max:	95 °C
Internal Min:	5 °C
External Max:	95 °C
External Min:	5 °C
Fluid:	H ₂ O (Water)

Fig. 21

Then the working screen appears and the device is ready for operation.

A	
0.00	
int °C	
- - -	- - -
ext °C	Safe Temp °C
int	5.00
	set °C
	2000
	Pump rpm

Fig. 22

Change the temperature setting with left knob (P).

Change the pump speed setting with the right knob (R).

In standby status, activate the tempering function by pressing the left knob (P). The device start tempering function. Meanwhile, the pump starts to run at once.

After pressing the power switch (11) the start-up of the compressor is delayed (about 3 minutes).

In working status, press the right knob (R) to stop operation of the pump. The tempering function and the pump stops.

Note: In standby status, press the right knob (R) to start the pump function, the tempering function will not be activated.

In working status, press the left knob (P) to stop the tempering, the pump keep running.

When the **WiCo** is switched on without being connected to the **station**, the green LED bar (C) and the green Bluetooth® LED (G) on **station** lights up.

The control elements of the **WiCo** allow to be locked by pressing the lock button (V) about 2 seconds, so no accidental changes during operation are possible (lock symbol  appears in the display).

By pressing lock button (V) about 2 seconds again, the control elements are released (lock symbol  disappears from the display).

Note: In an emergency, the device function can be turned off by pressing the “safe STOP” key (E) at the front of **station**. The LED bar (C) changed into red color.

A message appears in the display indicates that the **station** was forced off. To restart, press the power switch (6) or On/Off button (A) off and on again.

If the Bluetooth® function of the **WiCo** is activated, the Bluetooth® icon  appears on the screen and the user can search the **WiCo** by pressing the searching button (F). Then, a beep is heard.

Useful information

The **station** is controlled via a **WiCo**. If the **WiCo** is mounted on the **station**, data is exchanged between the **station** and **WiCo** via the contacts (D). The screen of the **WiCo** displays the home icon .

The **WiCo** is equipped with a USB socket (Universal Serial Bus) with which the **WiCo** can be connected with **station**, the USB icon  appears on the screen.

If the **WiCo** is not connected with **station** via a USB cable, the data exchange between the **station** and the **WiCo** via Bluetooth®. In this case, the Bluetooth® icon  is displayed.

Depending on the structure of the building, the **WiCo** can be operated at a distance up to 15 m from the **station**, using the Bluetooth® connection.

The **WiCo** could be either installed on the **station** or put on a safe place where is accessed easily by the user during operation.

If the **WiCo** is mounted on the **station**, the battery is charged through the contacts (D).

The battery could also be charged via the USB port on the **WiCo** (See “Charging the Battery Pack RB 1” in “Preparations” section).

Working screen at the time of delivery:




Fig. 23

Note: The wireless symbol  appears only when the **station** is switched on.

Explanation of symbols on the working screen:

The symbols displayed change depending on the status and settings of the **WiCo** and the **station**. The screen below shows the most significant symbols on the working screen.




Fig. 24

Bluetooth®:

This symbol means the **station** and the **WiCo** are communicating via Bluetooth®. The symbol no longer appears if no Bluetooth® communication is being performed.

99 Device number:

This figure appears when the "Device Number" is activated in menu option "Display". The figure can be set from 0 to 99. See section "Display" under "Menu (Details)".

Lock:

This symbol means that the function of the buttons and the rotary knobs for controlling the **WiCo** are disabled. The symbol no longer appears if the functions are enabled once again.

Temperature Sensor:

This symbol appears when the external temperature sensor is connected.

A Operating Mode:

This symbol indicates the operating mode currently selected (A, B, C, D).

USB:

This symbol means the **WiCo** is communicating or the battery is charged via a USB cable. The symbol no longer appears if no USB cable is being used for communicating with the **station**.

Home:

This symbol means that the **WiCo** is connected to the **station** and is communicating with the **station** via the contacts. The symbol no longer appears if the **WiCo** is removed from the **station**.

Battery pack:

This symbol indicates the charging status of the RB 1 battery pack within the **WiCo**.

The charging symbol appears if the **WiCo**

- is connected to a PC via a USB cable
- is connected to a **station** via a USB cable
- is connected to the power supply unit OS 1.0 via a USB cable
- is connected to the **station** via the charger contacts.

Heating

This symbol indicates that the heating function is active.  →  →  indicate active heating process.

Refrigerating

This symbol indicates that the refrigerating function is active.  →  indicates active refrigerating process.

Fluid level

This symbol indicates fluid level.

The red symbol  means the fluid is above the maximum fluid level. Excess fluid should be drained out. The red symbol  means the fluid is below the minimum fluid level. Please add fluid.

Pump:

This symbol indicates that the pump is activated.

Warning:

This symbol indicates that warning is active.

PC PC Control:


This symbol means that either the **station** or the **WiCo** is connected to a computer and is being controlled from the computer.

PR Program Control:

This symbol indicates that the **WiCo** is controlled by a program (see "PROGRAMS").

Menu navigation and structure:


Menu navigation:



- Press the "Menu" button (U).
- Select the menu by turning the right rotating/pressing knob (R) to the right or left.
- Open the menu item by pressing the right rotating/pressing knob (R).
- Turn the rotating/pressing knob (R) to select the desired menu option and edit the values or settings.
- Press the rotating/pressing knob (R) to get into sub menu items to active/inactive to switch settings or to confirm settings ("OK").
- Press the "Back" button (T) to a setting or cancel to return to the previous menu.
- Press the "Menu" button (U) to return directly to the working screen.

Note: When the tempering or pump functions are active, the menu is locked. On the display, the active menu option is highlighted in yellow. The active status of a menu item is marked with a check (✓).

Menu navigation: → Press the "Menu" button (U) and rotate rotating/pressing knob (R).
 ← Press the "Back" button (T) or the "Menu" button (U).



Note: If you press the "Menu" button (U), the system skips directly back to the working screen.
 If you press the "Back" button (T), the system skips back to the previous screen.

Fig. 26

Menu structure:

		Factory Settings	
Menu	TEMPERING	Control	
		Control Mode	Internal (int) activated
			External (ext) -
		Control Parameters	
		Automatic	Accurate activated
			Fast -
		Manual	Internal (Kp, Ti, Td, Ts, Prop_Bp, Prop_Bn) 60.0, 5.0, 0.3, 3s, +1.00K, -1.00K
			External (Kp, Ti, Td, Ts, Prop_Bp, Prop_Bn) 1.0, 15.0, 0.0, 90s, +1.00K, -1.00K
			Information (Kp, Ti, Td, Ts, Prop_Bp, Prop_Bn) yes
		Fluids	CF.EG28.N10.80.8 -
			CF.EG39.N20.80.16 -
			CF.EG44.N25.80.19 -
			CF.EG48.N30.80.22 -
			UF.Si.N30.150.10LV -
			HF.Si.20.200.50 -
			HF.Si.20.250.50A -
			H ₂ O (Water) activated
			Customized 1 -
			Customized 2 -
	Limits	Cooling	 activated
		External (ext)	Maximum 95 °C
			Minimum 5 °C
		Δ T (int - ext)	 100 °C
	Temperature Sensor	Power Output	Heating 100 %
		Calibration	Internal 2 Point Calibration
			External 3 Point Calibration
		Reset Calibration	Internal 2 Point Calibration
			External 3 Point Calibration
			 -
	Notification	ON	 -
		Hysteresis	 ± 0.0 K
		Sound	Infinite -
			Beep activated
PUMP	Limits	Speed	Maximum 4000 rpm
			Minimum 2000 rpm
		Pressure	Maximum 0.5 bar
TIMER	Set	Display	00:00 [hh:mm]
			 -
MODE	A		 activated
			 -
			 -
			 -

	DISPLAY	Main Screen	Pump	activated
			Timer	0
		Device number		
	GRAPH	Automatic		activated
				activated
		Manual	Axis Assignment	activated
			Set Temperature	activated
		Axis Scaling	Actual Internal Temp.	activated
			Actual External Temp.	activated
			X - Axis	15 min
				30 min
				60 min
				120 min
			Y - Axis	Maximum
				Minimum
	PROGRAMS	Program 1		
		...		
		Program 10		
	MODULES	M1 - Valve	ON	
			Default Status	Open
				Closed
			Control method	Level
				Start
			Pump OFF = Valve OFF	Stop
				Level 2
			Temperature	Level 4
				Higher
				Lower
			Int > Set Temperature	95 °C
				5 °C
		M2 - Output	Alarm	Error
				Warning
				Default Status
			Open	Open
				Closed
			Switch	ON
				Default Status
			Open	Open
				Closed
			Control method	Level
				Start
			Pump OFF = Switch OFF	Stop
				Level 2
			Temperature	Level 4
				Higher
				Lower
			Int > Set Temperature	95 °C
				5 °C
		M3 - Input	Ext. Standby	ON
				Default Status
			Open	Open
				Closed
	OUTGAS	Speed		2000 rpm
		Interval		10 s
	SAFETY	Time out	Set	01:00 [mm:ss]
				activated
			Last set Value	
			Tempering	Set
		Password	OFF	activated
				Last set Value
			Last Actual Value	
				Set
	SETTING	Languages	English	activated
			Deutsch	
			...	
		Display	Background	Black
				White
			Brightness	Standard Mode
				Battery Mode
			Firmware Update Info	50 %
				40 %
		Sound	Key Tone	
		Factory Settings		
		Communication	Device name	Circulator WiCo
			Bluetooth	activated
		Information	Tempering	yes
			Pump	yes
			Safety	yes
			Software	yes
			Intro	yes
			Firmware Update Info	yes

Menu (Details):



TEMPERING:

1. Control:

1) Control Mode:

Internal (int):

The temperature is regulated according to the internal temperature sensor.

Extern (ext):

The temperature is regulated according to the external temperature sensor.

2) Control parameters:

Automatic:

The optimal control parameters for PID temperature control are determined automatically. This is the recommended mode. Selecting "Automatic" allows you to set the dynamics for temperature control.

Accurate: precise tempering without overshooting.

Fast: fast tempering with minimal overshooting.




Fig. 27 (heating curve in "Automatic" mode)

Manual:

The control parameters for PID temperature control can be set manually.

„Manual“ should only be used in the event of special temperature control requirements.

When "Manual" is selected, the following parameters can be set for "Internal (int)" and "External (ext)" temperature control:

Kp: Proportional coefficient

The proportional coefficient **Kp** is the controller amplification and determines how strongly the control deviation (the difference between the target temperature and actual temperature) directly affects the control variable (on-time of the cooling). **Kp**-values that are too large can lead to the controller overshooting.

Ti: Integral time

The integral time **Ti** (s) is the correction time and determines how strongly the duration of the control deviation affects the control variable. **Ti** compensates for an existing control deviation. A high **Ti** means a smaller and slower effect on the control variable. **Ti**-values that are too small can lead to instability of the controller.

Td: Differential time

The differential time **Td** (s) is the derivative time and determines how strongly the rate of change of the control deviation affects the control variable. **Td** compensates for rapid control deviations. A high **Td** means a smaller and slower effect on the control variable. **Td**-values that are too large can lead to instability of the controller.

Ts: Sampling time

The sampling time **Ts** (s) is the time interval over which the control deviation is determined and the respective control variable (dependent on **Kp**, **Ti** and **Td**) is calculated.

Ts must be adjusted to match the response characteristic (total of all time constants) of the closed loop controlled system, so that the control variable can deliver a uniform and measurable change in the control deviation. **Ts** values that are too small or too large can lead to instability of the controller.

Prop Bp: Proportional Band Maximum.

Prop Bn: Proportional Band Minimum.

The Proportional Band is the range below (Prop_Bp) and above (Prop_Bn) the set value in which the control output value is calculated via the difference between the actual and the set value and the PID parameters.

Examples of non-optimal settings:




Fig. 28

3) Fluids:

Under the option "Fluids", a variety of fluids can be selected. The selected fluid limits the setting range of the target temperature. See table in the section "Fluid (Standard information for IKA fluid)".

The maximum and minimum temperature values of the selected fluid can be set within these limitations.

4) Cooling:

This menu option allows you to activate/deactivate the cooling function.

2. Limits:

1) External (ext):

Under the option "External (ext)", the maximum and minimum temperature for external temperature control can be limited.

2) ΔT (int - ext) I:

Under " ΔT (int - ext) I", the maximum difference between the internal and external temperature can be set.

The function limits the absolute temperature difference between the external system and internal bath during the heating up and cooling down process. This can protect sensitive device like the glass reactors from thermal shock (e.g. glass break).

3) Power Output:

Under "Power Output", the maximum heating output can be set as a percentage of nominal heat output.

3. Temperature sensor:

1) Calibration:

The internal and the external temperature measurement can be calibrated and adjusted.

You can select 2-point calibration or 3-point calibration for internal and external measurement.

Calibration proceeding (example: 2-point calibration):




Fig. 29

2) Reset Calibration:

By "Reset Calibration", the calibration value for the internal or the external temperature sensor will be deleted.

4. Notification:

1) ON:

This menu option allows you to activate/deactivate the "Notification" function. A green check mark and a beep informing when the target value (hysteresis) has been reached.

2) Hysteresis:

This menu option allows you to set the hysteresis from 0 to ± 2.0 K.

3) Sound:

Infinite: Beep until you press the "Back" button.

Beep: Single beep.



PUMP:

Limits:

In "Limits" menu option, you are allowed to set the maximum and minimum speed to the pump, and also the maximum pressure.



TIMER:

1. Set:

You can set a target time (duration) in this menu option. When device functions are started normally, this time is displayed on the working screen. The device functions stop automatically once this time has expired. The running time is then displayed again on the display.

Note: To deactivate the target time specification, set the target time to 00:00.

2. Display:

Activate the timer display on the main display (working screen).



MODE:

1. Operating Mode A:

After power-on/power failure no automatic restart of functions.

2. Operating Mode B:

After power-on/power failure automatic restart of functions, depending on previous settings.

3. Operating mode C:

Set values (set in A or B) cannot be changed.

After power-on/power failure automatic restart of functions, depending on previous settings.

4. Operating mode D:

Confirmation request for set value changes, if functions are active. After power-on/power failure no automatic restart of functions.



DISPLAY:

1. Main Screen:

In the "Main Screen" menu option, you can specify what information will be displayed on the screen.

2. Device Number:

In the "Device number" menu option, you are allowed to edit the device number from 0 to 99 that will be shown on the main screen.

When using multiple units, it may be helpful to identify **WiCo** and **station**.

Note: Labeling device number to the **station**.



GRAPH:

In this menu, you can set the options for the time-temperature diagram.

1. Automatic:

The scaling of the temperature axis (Y-axis) is automatically determined, depending on the target temperature and the actual internal and external temperature.

The time axis (X-axis) is permanently scaled to 30 minutes.

2. Manual:

1) Axis Assignment:

The temperature values to be displayed can be selected.

2) Axis scaling:

Scaling of the time (X) and temperature axis (Y) can be selected or set.



PROGRAMS:

Under programs, 10 user-defined temperature-time profiles can be created. A program can consist of up to 10 segments. Once a program has been selected, the following options are available:

1. Start:

Starts the program upon request of loop mode.

1) Infinite loop:

Upon completion of the last segment, the program continues with the first segment until you end the program by stopping a device function.

2) Loop Count:

Indicates the total number of loop cycles until program end.

When you start the program, the WiCo remind you to confirm the starting of the PR (Program) control.

If you need to pause the PR control, press menu button.

Note: At the end of the program all device functions are switched off.

2. Edit:

Edit/change the selected yellow highlighted program.

1) Seg No.:

Segment number.

2) Ctrl. Sensor (int/ext):

Determines whether control is through the internal (int) or external (ext) temperature sensor.

3) Temp.:

Target temperature.

4) Ctrl.Mode (Time / +/- x.x K):

In "Ctrl.Mode Time" the target values and settings of the segment are valid for the duration indicated in the column "Time hh:mm".

Afterwards, the next program segment is automatically executed.

In "Ctrl.Mode +/- x.xx K", the hysteresis (tolerance) of the actual temperature to the target temperature is set (e.g. +/- 0.1 K). The target values and settings of the segment are valid until the actual temperature reaches the target temperature +/- hysteresis for the first time. Afterwards, the next program segment is automatically executed.

5) Pump rpm:

Target speed of pump.

6) M1 (ON/OFF):

MODULES M1-Valve:

OFF: M1 valve in initial state

ON: M1 valve in inverted initial state.

Note: The initial state of the M1 valve is defined in "MODULES M1-Valve" in "Default Status" as "Open" or "Closed".

7) M2 (ON/OFF):

MODULES M2-Output switch

OFF: M2-output switch in initial state.

ON: M2-output switch in inverted initial state.

Note: The initial state of the M2-Output switch is defined in "MODULES M2-Output" in "Default Status" as "Open" or "Closed".

8) Cool:

Activate/deactivate the cooling function.

OFF: Cooling function is deactivated.

ON: Cooling function is activated.

Edit: Edit/change the selected program parameters.

Delete: Delete the selected program segment.

Insert: Inserts a new program segment after the selected segment.

Save: Saves changes.

3. Delete:

Deletes the selected yellow background program.

OK: Confirm the process.

Cancel: Cancels the process.

4. View:

Temperature-time overview displays for the program with segments of the selected program.

Note: If hysteresis is set as "Ctrl.Mode +/-x.xx K" for one or more segments in the program, the duration of the program cannot be determined.

Press and turn the knob (R) to display the segment details.

Once the program has been started, the program no., segment no. (active/total) and the remaining duration of the segment or hysteresis are displayed in the graph.



MODULES:

In "Modules" the multifunction port (8) outputs and inputs can be configured.

The outputs M1 and M2 can be controlled via "PROGRAMS".

1. M1-Valve:

1) ON:

Switches the external valve to the active state (inverted initial state).

Note: In a started „Program“ the M1 segment settings have higher priority.

2) Default Status:

Define the default state (OFF) of the external valve as "Open" or "Closed". It depends on the valve type (normally open or closed).

3) Control method:

Level: In this menu option, you can set and activate the start/stop fluid level for the valve.

When the menu option is activated, the valve opens when the filling level reaches the set start level and closes when the filling level reaches set stop level.

Pump OFF = Valve OFF: When the menu option is activated, the valve turns off when the pump is turned off.

Temperature: Temperature-dependent control of the M1-valve.

By entering the boundaries of "Higher" and "Lower", the temperature range is set in which the M1-valve can be opened.

If you activate a rule "Int/ext> Set Temperature" is an automatic control of the M1-valve dependence on the target temperature.

By the "Hysteresis" value (set temperature hysteresis) the control stability can be optimized.

Note: Int/ext stands according to the selected control mode (internal (int) or external (ext)) for the actual temperature to be controlled.

2. M2-Output:

1) Alarm:

Activate the alarm output (switch contact).

Error: In case of an error the "M2-Output — Switch" is activated (ON, inverted initial state).

Warning: In case of a warning the "M2-Output — Switch" is activated (ON, inverted initial state)

When "Warning" is activated, "Error" is simultaneously activated (see section "Error codes").

Default Status: Defines the default status (OFF) of the "M2-Output — Alarm" as "Open" or "Closed". It depends on the alarm type (normally open or closed).

2) Switch:

Activate the switch output.

ON:

Switches output into the active state (inverted initial state).

Note: When the "Program" is started, the M1 segment settings have higher priority.

Default Status:

Defines the initial state (OFF) of the switch output as "Open" or "Closed". It depends on the switch type (normally open or closed).

Control method:

Level: In this menu option, you can set and activate the start/stop fluid level for the switch.

When the menu option is activated, the switch starts to work when the filling level reaches the set start level and stops to work when the filling level reaches set stop level.

Pump OFF = Switch OFF: When the menu option is activated, the switch turns off when the pump is turned off.

Temperature: Temperature-dependent control of the M2-switch. By entering the boundaries of "Higher" and "Low", the temperature range is set in which the M2-switch can be opened.

If you activate a rule "Int/ext> Set Temperature" is an automatic control of the M2-switch dependence on the target temperature.

By the "Hysteresis" value (set temperature hysteresis) the control stability can be optimized.

Note: Int/ext stands according to the selected control mode (internal (int) or external (ext)) for the actual temperature to be controlled.

3. M3-Input:

1) Ext. Standby:

External standby input to stop the device functions Temper and Pump.

ON: Activate the ext. standby function. Device functions are stopped in the event of an inverted initial state (ON) at the input.

Default Status: Define the default status (OFF) of the input as "Open" (high level) or "Closed" (low level).



OUTGAS

In the menu option, you are allowed to set the outgas speed from 2000 rpm to 4000 rpm and outgas interval from 10 seconds to 240 seconds.

This function can be used when filling external devices such as laboratory reactors.



SAFETY

1. Time out:

1) Set:

In the menu option "Set", you can determine a time limit in the event of a communication breakdown between the **station** and the **WiCo**. The **station** continues to work with the preset target values until the preset time value has expired. Following that, the **station** runs using the preset safety temperature and safety speed.

Note: The initial time out is 30 seconds and you can define up to 60 minutes for this time limit.

2) Pump:

OFF: When this menu option is activated, the pump function switches off in the event of a communication breakdown between the **station** and the **WiCo**.

Last Set Value: When this menu option is activated, the pump continues to run with last set speed value in the event of a communication breakdown between the **station** and the **WiCo**.

Set: When this menu option is set and activated, the pump continues to run with the set speed value in the event of a communication breakdown between the **station** and the **WiCo**.

3) Tempering:

OFF: When this menu option is activated, the tempering function switches off in the event of a communication breakdown between the **station** and the **WiCo**.

Last Set Value: When this menu option is activated, the tempering function continues to work with last set temperature value in the event of a communication breakdown between the **station** and the **WiCo**.

Last Actual Value: When this menu option is activated, the tempering function continues to work with last actual temperature value in the event of a communication breakdown between the **station** and the **WiCo**.

Set: When this menu option is set and activated, the tempering function continues to work with the set temperature value in the event of a communication breakdown between the **station** and the **WiCo**.

2. Password:

In the "Password" menu, the menu settings can be locked by a 3-digit password.

3. Safe Temp. Confirmation:

This menu allows you to confirm the safe temperature of the **station** at start-up. This menu can confirm (safety temperature of the **station**) at start-up of the set "Safe Temp". A check mark (✓) indicates that the function is activated.



SETTING

1. Languages:

The "Language" menu option allows you to select the desired language.

2. Display:

The "Display" menu option allows you to change the background color and brightness of the working screen.

You can also activate the "Firmware Update Info" screen in this menu.

3. Sound:

The "Sound" menu option allows you to activate/deactivate the key tone.

4. Factory Settings:

Select the "Factory settings" menu option by turning and pressing the rotary/push knob. The system requests confirmation to restore the factory settings. Pressing the "OK" button resets all the system settings to the original standard values set at dispatch from the factory (see "Menu structure").

5. Communication

1) Device Name:

The "Device Name" menu option allows you to edit the device name. This can relatively identify the paired **station**.

2) Bluetooth:

The "Bluetooth®" option allows the user to activate/deactivate the "Bluetooth®" function. A tick shows that the option is activated.

6. Information:

The "Information" menu option offers you an overview of the most important system settings of the device.

Interface and output

The device can be connected to a PC and operated with the laboratory software Labworldsoft® through the RS 232 port, USB port or USB port on the **WiCo**.

The device software can also be updated with a PC via the RS 232 port / USB port.

Note: Please observe the system requirements together with the operating instructions and help section included with the software.

USB interface:

The Universal Serial Bus (USB) is a serial bus for connecting the device to the PC. Equipped with USB devices can be connected to a PC during operation (hot plugging). Connected devices and their properties are automatically recognized.

USB device drivers:

Connect the IKA device through the USB data cable to the PC. The data communication is via a virtual COM port. Configuration, command syntax and commands of the virtual COM ports are as described in RS 232 interface.

From Windows 10 and onwards the standard Windows USB driver is automatically loaded and a COM port number is assigned (find details in Windows Device Manager: "USB Serial Port (COMxx)"). If you have problems with USB communication, first ask your IT system administrator whether access to the USB interface is restricted for data security reasons.

Device software update:

For device software update, visit **IKA** website **www.ika.com** and enter the "Service" menu. Download and run the Firmware Update Tool.

Find and click the Firmware Update Tool in your PC after installation. Register your E-mail and password.

Connect the device to your PC via the USB cable. Then, you can update the device software according to instructions of the Firmware Update Tool.

RS 232 interface:

Configuration

- The functions of the interface connections between the stirrer machine and the automation system are chosen from the signals specified in EIA standard RS 232 in accordance with DIN 66 020 Part 1.
- For the electrical characteristics of the interface and the allocation of signal status, standard RS 232 applies in accordance with DIN 66 259 Part 1.
- Transmission procedure: asynchronous character transmission in start-stop mode.
- Type of transmission: full duplex.
- Character format: character representation in accordance with data format in DIN 66 022 for start-stop mode. 1 start bit; 7 character bits; 1 parity bit (even); 1 stop bit.
- Transmission speed: 9600 bit/s.
- Data flow control: none
- Access procedure: data transfer from the stirrer machine to the computer takes place only at the computer's request.

The following applies to the command set:

- Commands are generally sent from the computer (Master) to the device (Slave).
- The device sends only at the computer's request. Even fault indications cannot be sent spontaneously from the device to the computer (automation system).
- Commands are transmitted in capital letters.
- Commands and parameters including successive parameters are separated by at least one space (Code: hex 0x20).
- Each individual command (incl. parameters and data) and each response are terminated with CR LF (Code: hex 0x0d hex 0x0A) and have a maximum length of 80 characters.
- The decimal separator in a number is a dot (Code: hex 0x2E).

The NAMUR commands and the additional specific **IKA** commands serve only as low level commands for communication between the device and the PC. With a suitable terminal or communications programme these commands can be transmitted directly to the circulator equipment. The **IKA** software package, Labworldsoft, provides a convenient tool for controlling circulating equipment and collecting data under MS Windows, and includes graphical entry features, for pump motor speed ramps for example.

NAMUR Commands	Function
IN_PV_1	Read the external actual temperature
IN_PV_2	Read the internal actual temperature
IN_PV_3	Read the safety actual temperature
IN_PV_4	Read the pump actual speed
IN_SP_1	Read the internal setting temperature (if 0: internal control) Read the external setting temperature (if 1: external control)
IN_SP_3	Read the safety setting temperature
IN_SP_4	Read the pump setting speed
IN_TMODE	Read temperature control 0: internal control 1: external control
OUT_SP_1 xxx	Set the internal setting temperature XXX (if 0: internal control) Set the external setting temperature XXX (if 1: external control)
OUT_SP_12@n	Set the WD safety temperature with echo of the set (defined) value.
OUT_SP_4 xxx	Set the pump speed XXX
OUT_SP_42@n	Set the WD-safety speed with echo of the set (defined) value.
OUT_TMODE_0	Set to Internal temperature control
OUT_TMODE_1	Set to External temperature control
OUT_WD1@n	Start the watchdog mode 1 and set the time for the watchdog to n (20...1500) seconds. Echo of the Watchdog time. During a WD1-event, the tempering and pump functions are switched off. This command needs to be send within the watchdog time.
OUT_WD2@n	Start the watchdog mode 2 and set the watchdog time to n (20...1500) seconds. Echo of the watchdog time. During a WD2-event, the set temperature is changed to the WD safety temperature and the pump set speed is set to the WD safety speed. This command needs to be send within the watchdog time.
RESET	Reset the PC control and stop the device functions.
START_1	Start the tempering function.
START_4	Start the pump function
STOP_1	Stop the tempering function
STOP_4	Stop the pump function

Connections between device and external devices:

PC 1.1 Cable:

This cable is required to connect RS 232 port (5) to a PC.





Fig. 30

USB 2.0 Cable A - micro B:


This cable is required to connect USB port to a PC.



USB micro B
Fig. 31

USB 2.0 cable micro A - micro B:

This cable is required to connect **WiCo** to **station**.



USB micro B
Fig. 32

Connection **WiCo** to **station**:




Fig. 33



Connection the device to PC:






Fig. 34

Multifunction port:




Fig. 35

Maintenance and cleaning

To avoid contamination, it is essential to check and change the bath fluid regularly.


If water is used for bath fluid, we recommend to add Water bath protective media. The protective media stops the growth of algae, bacteria and other microorganisms.

To keep the full cooling performance, the dust filter of the cooler must be checked regularly and cleaned if necessary.

- Switch off the device and disconnect power cable.
- Open the front venting grid (3).
- Remove the bracket.
- Clean condenser dust filter with a vacuum cleaner or wash the filter in the water and dry it before assembly.

Note: Don't touch the condenser surface with hard parts.

Cleaning:



Disconnect main plug prior to cleaning!

Use only cleaning agents which have been approved by **IKA** to clean the device: water containing tenside / isopropyl alcohol.

- Wear protective gloves while cleaning the devices.
- Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- Do not allow moisture to get into the device when cleaning.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with **IKA** that this method does not damage the device.

Spare parts order:

When ordering spare parts, please give:

- machine type
- serial number, see type plate
- item and designation of the spare part, see **www.ika.com**
- software version.

Repair:

Please send in device for repair only after it has been cleaned and is free from any materials which may constitute a health hazard.

For repair, please request the "Decontamination Certificate" from **IKA**, or download printout of it from the **IKA** website **www.ika.com**.

If you require servicing, return the device in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Error codes

Any malfunctions during operation will be identified by an error message on the display.

Proceed as follows in such cases:

- Switch off device using the main switch at the back of the device
- Carry out corrective measures
- Restart device

Error code	Effect	Cause	Solution
Error 01	Pump off Heating/refrigerating off	No external sensor	- Check this sensor
Error 02	Pump off Heating/refrigerating off	Motor over current (rate current)	- Reduce pump motor speed - Use fluid with lower viscosity - Check if the pump impeller is blocked
Error 04	Pump off Heating/refrigerating off	Motor hall signal missing	- Reduce pump motor speed - Use fluid with lower viscosity - Check if the pump impeller is blocked
Error 05 	Pump off Heating/refrigerating off	Too high liquid level (when the tempering/pump function is activated)	- Check the fluid level and buoyage
Error 06 	Pump off Heating/refrigerating off	Too low liquid level (when the tempering/pump function is activated)	- Check the fluid level and buoyage
Error 09	Pump off Heating/refrigerating off	Device internal temperature is too high	- Check the ambient temperature and let the device cool down - Check the fan and clean the grids at the rear side
Error 11	Pump off Heating/refrigerating off	Temperature difference between control sensor and safety sensor is too much	- Check safety temperature circuit and bath fluid
Error 12	Pump off Heating/refrigerating off	Safety temperature alarm	- Check the bath temperature measurement
Error 13	Pump off Heating/refrigerating off	Heater switched off by safety circuit	- Check safety temperature set value, fluid level
Error 15	Pump off Heating/refrigerating off	Safety shutdown of the compressor	- Restart the device
Error 16	Pump off Heating/refrigerating off	Errors at the compressor (suction side)	- Restart the device
Error 17	Pump off Heating/refrigerating off	Errors at the compressor (pressure side)	- Restart the device
Error 20	Pump off Heating/refrigerating off	Fan cooling system error	- Restart the device

Warning message	Effect	Cause	Solution
High fluid level warning ⓘ	Warning message and high level icon appear	Too high fluid level	- Check the fluid level and buoyage - Drain excess fluid from bath.
Low fluid level warning ⓘ	Warning message and low level icon appear	Too low fluid level	- Check the fluid level and buoyage - Add more fluid to bath.

ⓘ Note: When the tempering or pump function is not activated, the screen show only these warning information if fluid level is too high or too low.

When the tempering or pump function is activated, if the fluid level goes slowly up/down to the warning level, the screen show high/low fluid level warning message and icon.

If the fluid level is corrected, the warning message and high/low level icon disappear.

If the fluid level increases/decreases further, the screen show Error 5/Error 6.

If the actions described fails to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department
- Send the device for repair, including a short description of the fault.

Accessories

Tubing and hoses:

LT 5.30	Insulated high temperature hose set, 800 mm, 2 pcs (M16x1)
LT 5.31	Insulated high temperature hose set, 1200 mm, 2 pcs (M16x1)
LT 5.32	Insulated high temperature hose set, 1600 mm, 2 pcs (M16x1)
H.PVC.8	PVC tube (nominal width 8)
H.PVC.12	PVC tube (nominal width 12)
H.SI.8	Silicone tube (nominal width 8)
H.SI.12	Silicone tube (nominal width 12)
H.PUR.8	PUR tube (nominal width 8 mm)
H.PUR.12	PUR tube (nominal width 12 mm)
H.FKM.8	FKM tube (nominal width 8 mm)
H.FKM.12	FKM tube (nominal width 12 mm)

Tubing insulations:

ISO. 8	Insulation (8 mm)
ISO.12	Insulation (12 mm)

Valve:

MV 1	Magnetic valve
CO V 1	Closed pressure valve

Additional accessories:

PC 1.1	Cable (RS 232)
Pt 100.30	Temperature sensor
Labworldsoft®	

See more accessories on www.ika.com.

Technical data

station		
Operating voltage	VAC	230 ± 10 % / 115 ± 10 % / 100 ± 10 %
Frequency	Hz	50/60
Max. input power	W	1800 (230 VAC) / 1500 (115 VAC) / 1210 (100 VAC)
Working temperature range	°C	- 30 ... + 100
Operating temperature range	°C	- 30 ... + 100
Temperature stability – internal temperature control, water (according to DIN 12876): -10 °C +70 °C	K	± 0.1 ± 0.05
Temperature control		PID (Automatic / user setting)
Connection for external temperature sensor		Pt 100
Temperature measurement, absolute accuracy internal (adjustable by calibration)	K	± 0.5
Temperature measurement, absolute accuracy external (adjustable by calibration)	K	± 0.5
External Pt 100.30 temperature sensor tolerance to DIN EN 60751 class A, ≤+/- (0.15+0.002x T), e.g. at max. 100°C (adjustable by calibration (ext))	K	± 0.35 (at 100°C)
Temperature setting		Knob on WiCo
Temperature setting resolution	K	0.1
Temperature display		TFT LCD on WiCo
Temperature display resolution	K	0.01
Classification according to DIN 12876-1		Class III (FL) suitable for flammable and non-flammable fluids
Safety circuit (adjustable)	°C	0 ... + 110
Safety temperature display		TFT LCD on WiCo
Heating capacity	W	1500 (230 VAC) / 1200 (115 VAC) / 910 (100 VAC)
Cooling capacity according to DIN 12876 (at 4000 rpm): + 20 °C + 10 °C 0 °C - 10 °C - 20 °C - 30 °C	W	620 600 480 350 200 80
Refrigerant		R 290 
Refrigerant quantity	g	79
Max. refrigerating system pressure	bar	21
Pump speed (adjustable)	rpm	2000 ... 4000
Max. pump pressure/suction	bar	0.5/0.25
Max. flow rate (at 0 bar)	l/min	21
Bath volume	l	1.4 ... 4.0
Fluid maximum viscosity	mm²/s	50
Low/high fluid level protection		Yes
Interface		USB, RS 232, multi-function port
Pump connection		M16x1
Permitted on-time	%	100
Noise level	dB (A)	58

IP code according to EN 60 529		IP 21
Protection class		I
Excess voltage category		II
Contamination level		2
Permitted ambient temperature	°C	+ 5 ... + 32
Permitted ambient humidity	%	80
Dimension (W x D x H)	mm	220 x 525 x 480
Weight	kg	26.5
Operation at a terrestrial altitude	m	max. 2000

WiCo

Permitted on time	%	100
Max. communication distance (dependent on the building)	m	15
Dimensions (W x D x H)	mm	160 x 40 x 105
Weight	kg	0.3
Ambient temperature	°C	+ 5 ... + 40
Ambient humidity (relative)	%	80
IP code according to EN 60529		IP 40
Interface		USB

RB 1 Battery pack

Voltage	V	3.7
Charging capacity	mAh	2000
Charging time	h	4.5
Working time	h	15
Battery type		Lithium-polymer

Note: Refrigerant must be disposed of in accordance with local and national regulations.

Subject to technical changes!

Warranty

In accordance with IKA Terms and Condition of sale, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Pump performance curve

Pump performance curve measured with water:

(Measurements done according DIN 12876-2 with water at 20 °C; pump in a closed-loop circuit).




Fig. 36



designed for scientists

IKA-Werke GmbH & Co. KG

Janke & Kunkel-Straße 10, 79219 Staufen, Germany

Phone: +49 7633 831-0, Fax: +49 7633 831-98

eMail: sales@ika.de

USA

IKA Works, Inc.

Phone: +1 910 452-7059

sales@ika.net

KOREA

IKA Korea Ltd.

Phone: +82 2 2136 6800

sales-lab@ika.kr

BRASIL

IKA Brasil

Phone: +55 19 3772-9600

sales@ika.net.br

MALAYSIA

IKA Works (Asia) Sdn Bhd

Phone: +60 3 6099-5666

sales.lab@ika.my

CHINA

IKA Works Guangzhou

Phone: +86 20 8222 6771

info@ika.cn

POLAND

IKA Poland Sp. z o.o.

Phone: +48 22 201 99 79

sales.poland@ika.com

JAPAN

IKA Japan K.K.

Phone: +81 6 6730 6781

info_japan@ika.ne.jp

INDIA

IKA India Private Limited

Phone: +91 80 26253 900

info@ika.in

UNITED KINGDOM

IKA England LTD.

Phone: +44 1865 986 162

sales.England@ika.com

VIETNAM

IKA Vietnam Company Limited

Phone: +84 28 38202142

sales.lab-vietnam@ika.com

THAILAND

IKA Works (Thailand) Co. Ltd.

Phone: +66 2059 4690

sales.lab-thailand@ika.com

TURKEY

IKA Turkey A.S

Phone: +90 216 394 43 43

sales.turkey@ika.com

KENYA

IKA Works Kenya Ltd.

Phone: +254 112 323 745

sales.kenya@ika.com

UGANDA

IKA Works Kampala Limited

Phone: +254 112 323 745

sales.uganda@ika.com

SPAIN

IKA Works Spain, S. L.

Barcelona

sales.spain@ika.com

Discover and order the fascinating products of IKA online:
www.ika.com



www.ika.com



IKAworlwide // #lookattheblue

2000050211_EN_HRC 2 GREEN control_062025_web